Instructions for Use

341712-G

English

OC-SENSOR DIANA Latex Reagent OC-SENSOR DIANA Buffer

REF V-PZ01
REF V-PZ03

INTENDED USE

For quantitative measurement of haemoglobin in faeces using automated immunochemical analysers [OC-SENSOR DIANA, OC-SENSOR PLEDIA].

INTRODUCTION

The amount of haemoglobin in faeces increases with illnesses that are accompanied by hemorrhagic lesions in the digestive tract, particularly in the lower digestive tract. Therefore measuring the amount of haemoglobin in faeces is an effective means of screening for early detection and treatment of colon cancer and other lower digestive tract illnesses that are accompanied by hemorrhaging¹⁾⁻⁸⁾.

Currently, the methods used for detection of haemoglobin include immunoassay latex agglutination⁹⁾, RPHA, EIA, and immunonephelometry. OC-SENSOR DIANA Latex Reagent is an immunoassay test reagent and used together with an automatic analyser for measurement of haemoglobin in faeces. OC-SENSOR DIANA Latex Reagent contains the latex particles sensitized to anti-human haemoglobin A₀ (HbA₀) polyclonal antibodies for optical measurement of the latex agglutination reaction.

PRINCIPLE OF THE METHOD

The test method is based on a latex agglutination reaction. A latex reagent is prepared by sensitizing anti-human HbA₀ antibodies to polystyrene latex particles. When this reagent is mixed with the sample, the anti-human HbA₀ antibodies, which were sensitized to latex, react with the haemoglobin in the sample, and the latex aggregate is formed in the latex agglutination reaction. The change in absorbance per unit time resulting from the latex agglutination reaction is proportional to the concentration of haemoglobin in the sample. A dose response curve of the absorbance unit (OD) vs. concentration is generated using the results obtained from the calibrators. The concentration of haemoglobin in the patient sample is determined from this curve.

REAGENTS

Reagents are stable until the date printed on the label assuming the container remains unopened at a storage temperature of 2-10°C.

OC-SENSOR DIANA Latex Reagent15mL, 5 vials, maximum 250 tests/vial
(1 vial contains 3mL of latex particles suspension coated with anti-human HbA₀ rabbit IgG)

OC-SENSOR DIANA Buffer500mL, 1 vial, maximum 1250 tests/vial
(1 vial contains 5.95g of N-2-Hydroxyethylpiperazine-N'-2-ethanesulfonic acid (HEPES))

Note: The above reagents are sold separately.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Before performing measurement, be sure to read the instruction manual of the analyser.
- Fully prepare and adjust the analyser before measurement.
- Create a calibration curve for each measurement. Also be sure to create a new calibration curve when OC-SENSOR DIANA Latex Reagent from a different lot is used.
- Use the Calibrator (REF V-PH51, V-PH52, V-IX50) for the creation of calibration curve.
- If the measurement results exceed the measurement range, use Diluent of Calibrator to dilute the sample and perform measurement again.
- Store all reagents at 2-10°C and avoid freezing.
- Do not use reagents that have passed their expiration date.
- The test sample may contain microorganisms. Therefore, use caution when handling. After use, all samples and other materials must be considered as medical waste and disposed of accordingly.
- Example of treatment: Soak for 1 hour or longer in a sodium hypochlorite solution (available chlorine concentration 1000 ppm or greater). (Neutralize any substances that contain acids before soaking.) Alternatively, treat in an autoclave at 121°C for 20 minutes. (Do not treat any items which sodium hypochlorite have adhered to in this way.)
- Dispose of used reagents and containers as medical waste in accordance with local regulations.
- If the product is used in any way other than that specified here, the reliability of measurement results cannot be guaranteed. Be sure to follow the procedure.
- A clinical diagnosis based on the measurement results must be a comprehensive judgment made by the attending physician, including factors such as clinical symptoms and other test results.

SAMPLE COLLECTION DEVICE

- Use the following designated sampling container for collecting specimens.
OC-Auto Sampling Bottle 3 (REF V-PZ25, V-PZ26)
- Collect faecal samples by scraping the surface of the stools in different areas. Collect the amount necessary to cover the groove of the probe.
- Check that the test sample has become fully suspended in the buffer solution inside the sampling bottle.
- Performance testing with the sample collection device demonstrated that samples stored at 2-10°C for 28 days have 95±14.7%, at 25°C for 7 days have 96±20.4%, for 14 days have 93±23.5%, and at 30°C for 7 days have 89±20.5%, for 14 days have 84±23.6% of haemoglobin recovery (in-house data, recovery rate shown as mean±2SD). However, the haemoglobin in some samples may undergo rapid denaturation or degradation, resulting in false negatives. Therefore, samples should be stored at 2-10°C and analysed as soon as possible.

PREPARATION OF REAGENTS

- OC-SENSOR DIANA Latex Reagent
Ready to use.
Before measurement, bring reagent to room temperature and invert the OC-SENSOR DIANA Latex Reagent gently several times to assure uniform suspension.
After opening, reagent stability (Decreasing DA values ≤5%):
● In refrigerator (2-10°C): 14 days
● Stored onboard : 7 days (Reagent cap needs to be closed when not testing)
Do not mix reagent with another bottle.
- OC-SENSOR DIANA Buffer
Ready to use.
Before use, stabilize the buffer to room temperature (20-25°C).
After opening, buffer stability:
● In refrigerator (2-10°C): 2 months.
● Stored onboard : 1 month.
Do not mix buffer with another lot.
Do not mix buffer from different bottles.
Do not use precipitated buffer, precipitate of buffer may cause problems.

REQUIRED MATERIALS, NOT PROVIDED BY THE MANUFACTURER (Prepare these materials before analysis)

- Wash solution: Sodium hypochlorite 0.15% (0.10%–0.30 % is acceptable)
- Purified water (for washing): Distilled or de-ionized water (1.0-10.0MΩ·cm is acceptable)
- Sample cup: To use sample cups, input the specification of the sample cup into the analyser.
- Printer Paper: Use thermal printer paper which fits the analyser.

TEST PROCEDURE

Perform measurements according to the instruction manual of the analyser.

- Set the OC-SENSOR DIANA Latex Reagent and OC-SENSOR DIANA Buffer on board the analyser.
- Input lot number parameters into the analyser.
- Check water, wash solution and drain tank volumes.
- Load control materials and/or sample collection devices.
- Start analysis.
- After measurement is completed, ensure the results of controls are within the acceptable range.

RESULTS

The reaction of each sample is compared to the calibration curve that was created previously. The concentration of HbA₀ (ng/mL) is found and a positive or negative judgment is made.

LIMITATIONS

- The test has not been validated for testing of patients with haemoglobinopathies.

ANALYTICAL SPECIFICITY

In comparison with Haemoglobin A₀, cross-reactivity of human haemoglobin Hb-S and Hb-C are 99% and 98%, respectively.

Cross-reactivity of 5 animal haemoglobin types such as bovine, swine, equine, rat and rabbit is approximately 2% or less, except for the equine, which shows 3.9%. However, no practical effect on the measurement is expected since it has been established that most haemoglobin contained in meat is digested and deactivated when it is consumed.

Almost no effect on the measurement value was found from the following co-existing substances at the indicated concentrations.

Substance	Concentration	Substance	Concentration
bilirubin	25mg/dL	Protein(bovine serum albumin)	2.5g/dL
lipids(intralipid)	0.6%	glucose	4.0g/dL
ascorbic acid	40mg/dL	barium sulfate	25mg/dL

INTERNAL QUALITY CONTROL

Each laboratory should establish a quality control program to monitor the performance of OC-SENSOR DIANA Latex Reagent. It is recommended to use the following materials for the quality control in your laboratory.

OC-Control LV1	(REF V-PH53)	OC-Control Low	(REF V-PH57)
OC-Control LV2	(REF V-PH54)	OC-Control High	(REF V-PH58)
OC-Control LV3	(REF V-PH59)		

PERFORMANCE CHARACTERISTICS

- Sensitivity
When standard solution of 50ng/mL HbA₀ was measured in comparison to the blind reagent test, the change in absorbance that was obtained was 0.002 or more.
- Accuracy
When control samples of known concentration were measured, the value obtained was within ±15% of the indicated value.
- Inter/Intra-day reproducibility
When the same sample was measured 10 times within the same run, the coefficient of variation (CV) for the values obtained was 5% or less.
When the same control and the same faeces samples (129ng/mL-630ng/mL) were measured for 20days, the coefficient of variation (CV) for the values obtained for inter/intra-day was within 7%.
- Measurement range
50 – 1000ng/mL (10µg/g – 200µg/g) (faeces)
- Acceptable range between-reagent batches
- OC-SENSOR DIANA Latex Reagent
When control samples with known concentration were measured, the value obtained was within 100±15% of the indicated value.
- OC-SENSOR DIANA Buffer
When OC-Calibrator 1 kit was measured, the OD value obtained was within 100±5% as compared with the OD value obtained by measuring OC-Calibrator 1 kit using in-company standard products.
- Cut-off value should be set on the demand of each laboratory¹⁰⁾.

PRODUCT CODE, PRODUCT NAME & STORAGE

Product code	Product name	Contents	Storage	Analyser
V-PZ01	OC-SENSOR DIANA Latex Reagent	5×15mL	2-10°C	DIANA, PLEDIA
V-PZ03	OC-SENSOR DIANA Buffer	1×500mL	2-10°C	DIANA, PLEDIA
V-PH51	OC-Calibrator 1 kit	1×3mL	2-8°C	DIANA
V-PH52	OC-Calibrator 2 kit	1×3mL	2-8°C	PLEDIA
V-IX50	OC-Standard	10×1mL	2-10°C	DIANA
V-PH53	OC-Control LV1	2×5mL	2-8°C	DIANA, PLEDIA
V-PH54	OC-Control LV2	2×5mL	2-8°C	DIANA, PLEDIA
V-PH59	OC-Control LV3	2×5mL	2-8°C	DIANA, PLEDIA
V-PH57	OC-Control Low	5×1mL	2-10°C	DIANA
V-PH58	OC-Control High	5×1mL	2-10°C	DIANA
V-PZ25	OC-Auto Sampling Bottle 3	100 bottles	1-30°C	DIANA, PLEDIA
V-PZ26	OC-Auto Sampling Bottle 3 without barcode	100 bottles	1-30°C	DIANA, PLEDIA

REFERENCE

- J.S.Mandel et al.: New England Journal of Medicine, 328: 1365-1371, 1993.
- S.J.Winawer et al.: Journal of National Cancer Institute, 85 (16): 1311-1318, 1993.
- N.Hiwatachi et al.: Japan Journal of National Cancer Research, 84: 1110-1112, 1993.
- H.Saito et al.: International Journal of Cancer, 61: 465-469, 1995.
- B. Levin et al.: Gastroenterology, 134(5): 1570-1595, 2008.
- S.Ciatto et al.: British Journal of Cancer, 96: 218-221, 2007.
- G.Grazzini et al.: European Journal of Cancer Prevention, 13: 16-26, 2006.
- E. Fenocchi et al.: European Journal of Cancer Prevention, 15: 384-390, 2006.
- T. Takeshita et al.: Journal of Coloproctology, 38: 780-783, 1985.
- LGM Rossum et al.: British Journal of Cancer, 101: 1274-1281, 2009.

Advena Ltd.

Tower Business Centre, 2nd Flr., Tower Street,
Swatar, BKR 4013 Malta

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EIKEN CHEMICAL CO., LTD.